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## SYNTHETIC FIBER MOLECULAR ORIENTATION TESTER

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## ■Feature

The molecules of synthetic fiber are arranged to axis-direction in the process.

The degree of arrangement is called orientation.

The orientation influences so much to the fiber's quality.

There are some methods to measure orientation. We applied the practical birefringence ( $\Delta N$ ) method. Moreover, our original step-compensator method realized efficient measurement and has following feature.

### 1. Measuring the $\Delta N$

Measuring of  $\Delta N$  is possible by simple construction. It is only added the polarizer, retardation plate and visual processor to the biological microscope.

### 2. Measuring on the monitor

Image of microscope is inputted by the CCD camera, and displayed to the monitor.

Measuring of  $\Delta N$  is possible by mouse operation on the monitor.

### 3. Automatic processing of the statistics

Measuring data of  $\Delta N$  is processed automatically and average, standard deviation, co-efficient of variation are printed.

### 4. Application of the visual processor

Visual processor can apply to analyzing the diameter, section (round, irregular, hollow) of the fiber, and other various use.

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
## Principle of Measuring Orientation

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### ● Relation between Orientation and $\Delta N$

Case where the molecules are oriented to axis-direction, axis-direction's refractive index ( $N_x$ ) and diameter-direction's refractive-index ( $N_y$ ) are different, on that occasion, the birefringence ( $\Delta N$ ) occurs.

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
### ● Measuring $\Delta N$ with Step-compensator method

If fiber has birefringence ( $\Delta N$ ) character, interference stripe (rainbow color) occurs.  
Case where interference is compensated by combining the retardation plate, black-stripe occurs. The value of  $\Delta N$  is calculated automatically.

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## Specifications

Type	DELTA-N
Measuring yarn	Synthetic fiber (round section)
Microscope	Biological microscope <ul style="list-style-type: none"> <li>● Objective (×4, ×10, ×20, ×40)</li> <li>● Eyepiece (×10)</li> <li>● Analyzer (with-in eyepiece tube)</li> </ul>
Retardation plate	1/4λ、1/2λ、1λ、3λ、5λ
Polarizer	Retardation plate setting type
CCD camera	Monochrome type : 1/3 inch Color type : 1/2 inch (Option)
Monitor	Monochrome type : 14 inch Color type : 14 inch (Option)
Visual processor	Soft-ware (△N, Length, Circle, Polygon, Free-curve, Angle, etc)
Printer	Thermal-type (40 letters/line)
Power	AC 100V (Standard), 50/60Hz Approx. 200W
Setting space	Approx. 1000(W)×600(D)×700(H) (mm)
weight	Approx. 40kg

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